



Test kit, chloride, model CDS-DT

Ürün Numarası:

2580600

TRY Fiyat:

Bize Ulaşın

3 gün içinde teslimat

Tehlikeli



The single-parameter Model CDS-DT uses the Hach Digital Titrator with Silver Nitrate Titrant Cartridges to determine chloride concentrations. Kit includes Digital Titrator, reagents and labware, manual and carrying case. Range: 10-10,000 mg/L. Contains approximately 100 tests.

Hach Digital Titrator offers higher levels of precision and accuracy than drop count procedures, with reproducibility comparable to a buret titration

Rugged titration device is ideal for field use, with no glass burets required

Teknik özellikler

Metot:	Digital Titrator/ Silver nitrate
Muhafaza Türü:	F
Ölçüm aralığı:	10 - 10,000 mg/L Cl^-
Ölçüm aralığı (2):	to
Parametre:	Chloride - As Cl^- , 10-10,000 mg/L Cl^-
Platform:	Dijital Titratör
Smallest increment:	10
Test kiti modeli:	CDS-DT
Test sayısı:	100

Gerekli Aksesuarlar

- Pipet, değişken hacimli, 1,0 - 5,0 mL (Item BBP065)
- Gloves, Disposable, Powder Free, Nitrile, Medium (Item 2550502)
- Gloves, Disposable, Powder Free, Nitrile, X-Large (Item 2550504)
- Gloves, Disposable, Powder Free, Nitrile, Small (Item 2550501)
- Gloves, Disposable, Powder Free, Nitrile, Large (Item 2550503)