



## Reagent and apparatus set, water in oil

Ürün Numarası:

2237300

Tehlikeli

**ÜRETİMDEN KALDIRILAN ÜRÜN**  
Bu ürün artık mevcut değil.



### Model WO-1 Test Kit

The Model WO-1 test kit provides critical data in oil quality determinations, oil line operation, and diesel engine maintenance. The kit has proven particularly useful on board ships. Reagent is packaged in unit-dose vials, and parts in contact with the sample are disposable, so clean-up is eliminated.

**Measures low levels of water in oil**

**Provides critical data in oil quality deterinations, oil line operation, and diesel engine maintenance**

**Reagent is packaged in unit-dose vials, and parts in contact with the sample are disposable, so clean-up is eliminated**

### Teknik özellikler

Metot:	Volumetric displacement/ Calcium hydride
Muhafaza Türü:	no case
Ölçüm aralığı:	0,05 - 1 %
Ölçüm aralığı (2):	0 - 10 %
Parametre:	Water-in-Oil - As %
Smallest increment:	0.05% 0.20%
Test kiti modeli:	WO-1
Test sayısı:	25