



## Test strips copper , 0-3 mg/L

Ürün Numarası:

2745125

TRY Fiyat:  
Mevcut

Bize Ulaşın

Tehlikeli



Copper exists in water as a soluble salt or as suspended solids. Concentrations exceeding 0.1 mg/L are useful for controlling algae and plankton growth. Quantities ranging from 0.02 to 0.1 mg/L are toxic for some fish, so its use for treating fish disease requires careful monitoring. With Hach Water Quality Test Strips, technicians in the field can test many samples in only a few minutes, and make immediate evaluations on-site. No measuring, set-up, clean-up, or chemical handling are necessary. Hach test strips are also used in laboratories all over the world for pre-test screening-to detect the presence of materials that might interfere with lab testing. Contains 25 tests.

**Easy to use, disposable, and inexpensive**

**A great way to obtain quick, quantitative answers in the field or in the lab**

### Teknik özellikler

Birim:	25 tests
Dipnot:	*mg/L unless otherwise noted; ppb = µg/L; ppm = mg/L.; gpg = grains per gallon; 1 gpg = 17.1 mg/L or 17.1 ppm.
En küçük artış adımları:	Steps: 0, 0.2, 0.5, 1, 3 ppm
Ölçüm aralığı:	0 - 3 ppm Cu
Parametre:	Copper
Platform:	Test Strip
Test kiti modeli:	Test Strips
Test sayısı:	25