



Membrane filters millipore S-pak type HC , 47 mm, 0.7 um pore size, pk/150

Ürün Numarası:

2936200

TRY Fiyat:

Bize Ulaşın

Mevcut

Certified for fecal coliform analysis. Larger pore size (0.7 um) and funnel-like pores creates a less taxing environment for stressed organisms: lessens media loss from evaporation and easier filtration of high particulate samples common to effluents. Whi

Teknik özellikler

Çap:	47 mm mm
Gözenek boyutu:	0.7 µm
Hacim / Paket boyutu:	150
Malzeme:	Mixed-Cellulose Ester