



Filter, acrocap, pore size 0.2 µm

Ürün Numarası: 2323810
TRY Fiyat: Bize Ulaşın
Mevcut

For sterilization or clarification of up to three liters of culture media or other aqueous solutions. Ideal for sterilization of solutions that tend to foam when filtered under vacuum. Inlet accepts 6.4 mm (0.25") tubing or most male luer fittings.

Teknik özellikler

Çap:	NA
Filter Type:	Filter, acrocap
Gözenek boyutu:	0.2 µm
Packing:	NA