




Filter, acrodisc, pore size 0.2 µm

Ürün Numarası:	2415600
TRY Fiyat:	Bize Ulaşın
Mevcut	

For filtration of dilute biological fluids. Sterilized by gamma irradiation to eliminate potential contamination by EtO residuals. Inlet/Outlet connections: 13/25 mm: Female Luer-Lok inlet, standard male slip luer outlet. Pore size: 0.2 µm. Diameter: 25mm.

Teknik özellikler

Çap:	NA
Filter Type:	Filter, acrodisc
Gözenek boyutu:	0.2 µm
Malzeme:	Polysulfone
Packing:	NA