



Microfil funnel & filter, 100 mL, 0.45 µm 47 mm white gridded (millipore)

Ürün Numarası:

2586300

TRY Fiyat:

Bize Ulaşın

Mevcut

Designed to eliminate cleaning and sanitization of filter funnels between uses. Graduated (20, 50, and 100 mL) 100 mL funnels are sterile and ready-to-use. Funnel creates leak-proof seal when seated onto Microfil Stand Alone Support. Excellent for testing water samples from potable, non-potable, source, ambient, industrial, food, beverage, pharmaceutical and wastewater sources. Microfil funnels come with 150 Millipore filters (white; gridded; cellulose and mixed ester composition; 0.45 µm pore size; 47 mm diameter; sterile; individually packed). Microfil V includes integrated filter for added convenience and can be used with standard No. 8 stopper on any filter holder or manifold.

Teknik özellikler

| | |
|-----------------------|-----------------------|
| Çap: | 47 mm |
| Gözenek boyutu: | 0.45 µm |
| Hacim / Paket boyutu: | 150 |
| Kapasite: | 100 mL |
| Malzeme: | Mixed-Cellulose Ester |