



Vacuum filter microfunnel, pk/50.

Ürün Numarası: 2831500
TRY Fiyat: Bize Ulaşın
Teslimat tarihi için arayın

Feature a completely assembled filtration unit: Funnel, filter, nutrient pad and petri dish. Filter: GN-6 Metrical with 0.45 µm pore size and 47 mm diameter; white, gridded surface. Funnel: Polypropylene Petri Dish: polystyrene lid and polypropylene base Nutrient Pad: Cellulose. Fits 2831500, 2846100 and 2846150 manifolds. Sterile and individually bagged.

Teknik özellikler

Çap: 47 mm
Gözenek boyutu: 0.45 µm
Hacim / Paket boyutu: 50
Malzeme: Mixed-Cellulose Ester