



Membrane filters millipore S-pak type HA, 47 mm, 0.45 um pore size, pk/600

Ürün Numarası: 2936500
TRY Fiyat: Bize Ulaşın
Mevcut

For membrane filtration in microbiological analysis of aqueous samples. "HA" signifies certified compliance with Standard Methods criteria with pore size of 0.45 um. White, cellulose and mixed-ester composition. S-Pak (filter only) is sterile. Diameter: 4

Teknik özellikler

Çap:	47 mm mm
Gözenek boyutu:	0.45 µm
Hacim / Paket boyutu:	600
Malzeme:	Mixed-Cellulose Ester