



Filter, membrane, black, pore size 0.45 µm, diameter 47 mm, 200/pk, w/out pads

Ürün Numarası:

2490300

TRY Fiyat:

Bize Ulaşın

2 hafta içinde teslimat

Provides contrast for microbiological and particulate analysis. Modified polyethersulfone creates dark background for excellent contrast useful in enumerating opaque colonies and light-covered particulate. White grid line aids enumeration. Certified for membrane filtration analysis. Sterile and individually packaged. Pore size 0.45 µm. Diameter: 47 mm.

Teknik özellikler

Çap:	47 mm
Filter Type:	Filter
Gözenek boyutu:	0.45 µm
Hacim / Paket boyutu:	200 /pk
Malzeme:	Modified polyethersulfone