



Filter, membrane, pore size 0.2 µm, diameter 47 mm, 200/pk (Gelman)

Ürün Numarası:

2618800

TRY Fiyat:

Bize Ulaşın

2 hafta içinde teslimat

For isolation and enumeration of stressed organisms (*Pseudomonas* species). High retention and flow rate. Certified for use in membrane filtration analysis. White polyethersulfone media. STerile and individually packaged. Pore size: 0.2 µm. Diameter: 47mm.

Teknik özellikler

Çap:	47 mm
Filter Type:	Filter
Gözenek boyutu:	0.2 µm
Hacim / Paket boyutu:	200 /pk
Malzeme:	Polyethersulfone