



## Sample vial, 120 mL w/ thiosulfate & custody tie, 200/pk

Ürün Numarası: 2634001  
TRY Fiyat: Bize Ulaşın  
Teslimat tarihi için arayın

### Sample Vial with Thiosulfate and Custody Tie

120 mL sterile, hinged-cap vials meet all EPA requirements for microbiology sampling. Collection, incubation, and accurate reading can all be done in the vial without a transfer step. 100 mL raised fill line. Polypropylene. ( **Please Note:** These sample vials are not autoclaveable)

### Teknik özellikler

|                       |                                 |
|-----------------------|---------------------------------|
| Birim:                | 200 /pk                         |
| Çap:                  | 43mm                            |
| Hacim / Paket boyutu: | 200 /pk                         |
| Kapasite:             | fill line at 100 mL (+/-2.5 mL) |
| Önemli özellik:       | Sample Vial, with Thiosulfate   |
| Uzunluk:              | 93mm                            |