



Agar LES, m-Endo, dehydrated, 500 g

Ürün Numarası:

2406634

Tehlikeli

TRY Fiyat:
Mevcut

Bize Ulaşın



Complete Solutions for Testing Recreational Water

Hach offers a variety of media and apparatus to meet all of your recreational testing needs.

Complete Solutions for Testing Wastewater

Hach has prepared media and convenient apparatus to make wastewater and non-potable water testing fast and easy.

Complete Solutions for Testing Drinking Water

Whichever method a customer chooses for examining drinking

water samples, Hach has prepared media and convenient apparatus to make

testing easy and fast.

Teknik özellikler

Hacim / Paket boyutu:

500 g

Media Type:

m-Endo LES Agar