



## Economy kit m-ColiBlue24, glass ampules, 200 test

Ürün Numarası: 2779200  
TRY Fiyat: Bize Ulaşın  
Teslimat tarihi için arayın

In 24 hours, Hach's m-ColiBlue24 Broth gives you simultaneous results for total coliforms and E. coli. Use m-ColiBlue24 Broth to accurately monitor and evaluate total coliform and E. coli in drinking water; wastewater; bottled water; beverages; surface, ground, and well water; and ultrapure, chemical processing and pharmaceutical processing waters.

With ready-to-use m-ColiBlue24 Broth, you can differentiate colonies easily - red and blue indicate total coliforms and blue specifies E. coli. Get optimal recovery of stressed and injured organisms. Minimize background growth of non-coliform bacteria.

m-ColiBlue24 should only be used with pads manufactured by the following companies/brands: Pall Corporation (e.g., Sterile Petri Dishes with Absorbant Pads), Private Label/OEM products manufactured by Pall Gelman, Sartorius (cellulose and glass fiber pads), Sartorius BioSart products with the integrated Petri dish, and Millipore EO sterilized pads. No other pad is compatible with the product, and the use of any pad other than those recommended may result in unreliable, inaccurate data.

**Read and confirm results in 24 hours**

**Enumerate total coliforms and E. coli on one petri dish**

**Conduct emergency testing and routine monitoring**

**Get superior sensitivity - 1 CFU/100mL**

**Eliminate special equipment and ultraviolet lamp**