



Ampules, m-TGE/TTC broth, plastic, pk/50

Ürün Numarası: 2428450
TRY Fiyat: Bize Ulaşın
Mevcut

m-TGE w/TTC, 2 mL plastic ampules. pk/50.

Complete Solutions for Testing Recreational Water

Hach offers a variety of media and apparatus to meet all of your recreational testing needs.

Complete Solutions for Testing Wastewater

Hach has prepared media and convenient apparatus to make wastewater and non-potable water testing fast and easy.

Complete Solutions for Testing Drinking Water

Whichever method a customer chooses for examining drinking

water samples, Hach has prepared media and convenient apparatus to make

testing easy and fast.

Teknik özellikler

Hacim / Paket boyutu:	50 /pk
Malzeme:	Plastic
Media Type:	m-TGE/TTC Broth