



Plastic ampules, pseudomonas broth, pk/50

Ürün Numarası:	2812250
TRY Fiyat:	Bize Ulaşın
Mevcut	

Pseudomonas broth plastic ampules, 2 mL. For detection of Pseudomonas species. pk/50.

Complete Solutions for Testing Recreational Water

Hach offers a variety of media and apparatus to meet all of your recreational testing needs.

Complete Solutions for Testing Wastewater

Hach has prepared media and convenient apparatus to make wastewater and non-potable water testing fast and easy.

Complete Solutions for Testing Drinking Water

Whichever method a customer chooses for examining drinking water samples, Hach has prepared media and convenient apparatus to make testing easy and fast.

Teknik özellikler

Hacim / Paket boyutu:	50 /pk
Malzeme:	Plastic
Media Type:	Pseudomonas Broth