



Plastic ampules R2A, pk/50

Ürün Numarası: 2812350
TRY Fiyat: Bize Ulaşın
Mevcut

R2A Broth Ampules; Plastic

R2A Plastic Ampules for detection of stressed heterotrophic bacteria; pk/50.

Complete Solutions for Testing Recreational Water

Hach offers a variety of media and apparatus to meet all of your recreational testing needs.

Complete Solutions for Testing Wastewater

Hach has prepared media and convenient apparatus to make wastewater and non-potable water testing fast and easy.

Complete Solutions for Testing Drinking Water

Whichever method a customer chooses for examining drinking water samples, Hach has prepared media and convenient apparatus to make testing easy and fast.

Teknik özellikler

Hacim / Paket boyutu:	50 /pk
Malzeme:	Plastic
Media Type:	R2A