




StablcalÂ Turbidity Standard, 0.5 NTU 500ml

Ürün Numarası:

2698049

Tehlikeli

TRY Fiyat:
Mevcut

Bize Ulaşın



Stablcal Turbidity Standard for all Turbidimeters. Stablcal Stabilized Formazin Turbidity Standards are true Formazin dilutions developed for use in any turbidimeter.

Precise, Pre-mixed Formazin Solutions

Stablcal Stabilized Formazin Turbidity Standards are true Formazin dilutions developed for use in any turbidimeter. With proprietary manufacturing technology, Hach prepares Stablcal Standards in precise concentrations as low as 0.060 NTU for conventional turbidimeters. Shelf life is guaranteed for two years from the date of manufacture. Whether you are monitoring wastewater effluent, high-quality treated drinking water or an industrial process stream, you can be completely confident in the results you obtain with Hach Stablcal Stabilized Formazin Standards.

No special glassware or difficult lab techniques

Hach was the first to apply Formazin as a turbidity standard in the drinking water industry, more than 40 years ago. Since then, Formazin has been the only primary standard that is accepted virtually throughout the world for regulatory reporting. All other turbidity standards, both alternate and secondary, must be traced to it.

Developed for use in any turbidimeter.

With proprietary manufacturing technology, Hach prepares Stablcal Standards in precise concentrations as low as 0.060 NTU for conventional turbidimeters.

Teknik özellikler

Hacim / Paket boyutu:	500 mL
Konsantrasyon:	0.5 NTU
Ölçüm aralığı:	N/A
Parametre:	Turbidity
Saklama koşulları:	Store between 5 - 25 °C, away from direct sunlight

Gerekli Aksesuarlar

- Gloves, Disposable, Powder Free, Nitrile, Medium (Item 2550502)
- Gloves, Disposable, Powder Free, Nitrile, X-Large (Item 2550504)
- Gloves, Disposable, Powder Free, Nitrile, Small (Item 2550501)