



9245/9240 Sodyum Standardı, 100 mg/L, 1 L

Ürün Numarası: 2834253
TRY Fiyat: Bize Ulaşın
Mevcut

Teknik özellikler

Cihaz:	9245 & 9240
Hacim / Paket boyutu:	1 L
Konsantrasyon:	100 mg/L
Parametre:	Sodium
Saklama koşulları:	Keep containers tightly closed in a dry, cool and well-ventilated place. Protect from moisture.

Gerekli Aksesuarlar

- Gloves, Disposable, Powder Free, Nitrile, Medium (Item 2550502)
- Gloves, Disposable, Powder Free, Nitrile, X-Large (Item 2550504)
- Gloves, Disposable, Powder Free, Nitrile, Small (Item 2550501)
- Gloves, Disposable, Powder Free, Nitrile, Large (Item 2550503)