



Recharge kit of 4 pre-filled cartridges with premounted 2956A membranes for A1100 oxygen sensors with 33051-xT cap only

Ürün Numarası:

2956A-AT05

Tehlikeli

TRY Fiyat:

Bize Ulaşın



Includes o-rings and cleaning tools. Special low concentration electrolyte.

Teknik özellikler

Cihaz:	Orbisphere A1100-S00T O2 sensor, Orbisphere GA2x00-S00T O2 sensor
Parametre:	O2
Saklama koşulları:	Store in a dry place at room temperature

Gerekli Aksesuarlar

- Gloves, Disposable, Powder Free, Nitrile, Medium (Item 2550502)
- Gloves, Disposable, Powder Free, Nitrile, X-Large (Item 2550504)
- Gloves, Disposable, Powder Free, Nitrile, Small (Item 2550501)
- Gloves, Disposable, Powder Free, Nitrile, Large (Item 2550503)