




Colour Standard Solution, 500 Pt Co Units, pk/16 - 10 mL Voluette Ampules

Ürün Numarası:

141410-EU

Tehlikeli

TRY Fiyat:
Mevcut

Bize Ulaşın



For accuracy checks in Color determinations. 500 platinum-cobalt units. pk/16 - 10 mL Voluette Ampules.

Simple Accuracy Checks

Voluette and PourRite Ampule Standards are designed to work directly with Standard Additions accuracy checks in Hach procedures and photometer software. Use the TenSette Pipet to add precise spikes of standard and check percent recovery, to confirm analytical system performance.

Minimize waste

Open a fresh standard ampule for each accuracy check. Each ampule contains enough standard for one set of standard additions.

Reliable results

Sealed ampules prevent accidental contamination or evaporation of standard solutions.

Teknik özellikler

Hacim / Paket boyutu:	16 x 10 mL ampules
Konsantrasyon:	500 Pt-Co units
Parametre:	Color
Platform:	Çözelti
Tanım:	Color Standards