



Be Right™



Test kit, nitrogen, nitrate color disc, model NI-14

Ürün Numarası:

1416100-EU

Tehlikeli

TRY Fiyat:

Bize Ulaşın

3 gün içinde teslimat



Utilizing the simple color disc method, this kit measures Nitrogen, Nitrate as NO₃ and contains 50 tests.

Simple, low-cost color disc method is the most accurate visual test kit method

Kit contains all required reagents and apparatus in a rugged carrying case for analysis in the field

Teknik özellikler

En küçük artış adımları:	0,05
Metot:	Color disc/ Cadmium reduction
Model:	NI-14
Muhafaza Türü:	B
Ölçüm aralığı:	0 - 1 mg/L NO ₃ ⁻ -N
Ölçüm aralığı (2):	0 - 10 mg/L NO ₃ ⁻ -N
Parameter:	Nitrogen, Nitrate
Parametre:	Nitrogen, Nitrate - As NO ₃ ⁻ -N
Platform:	Renk Diski
Smallest increment:	0,05
	0,5
Standart aralık:	0 - 1 mg/L NO ₃ ⁻ -N; 0 - 10 mg/L NO ₃ ⁻ -N
Test kiti modeli:	NI-14**
Test sayısı:	50