



Be Right™



NitraVer 5 Nitrate Reagent Powder Pillows, 5 mL, pk/100

Ürün Numarası:

1403599-EU

Tehlikeli

TRY Fiyat:
Mevcut

Bize Ulaşın



Powder reagent in sealed packages for the determination of Nitrogen-Nitrate up to 50 mg/l NO₃-N.

For high range nitrate determination by the NitraVer 5 cadmium reduction method. Range: up to 50 mg/L as nitrate-N. Sample Size: 5 mL. Pack of 100 powder pillows.

Teknik özellikler

Hacim / Paket boyutu:

100 /pk

Metot:

Cadmium Reduction

Platform:

Powder Pillows